

Penicillin-Streptomycin, liquid
Prepared with 10,000 units/mL Penicillin
G Sodium and 10,000 ug/mL Streptomycin

Catalog Number: 1931619

Lot Number: 2214193

Storage Temperature: -5 to -20C

Expiration Date: 2023-03-31

For research use only. Not for use in diagnostic procedures.

| TEST | SPECIFICATION | RESULT | UNITS |
|-------------------------|----------------|----------|---------|
| ¹ Osmolality | >=310 to <=350 | 324 | mOsm/kg |
| ² pH | >=4.0 to <=7.0 | 5.4 | |
| ³ Sterility | Negative | Negative | |

Read SDS.



Quality Systems Department

Date: 30-Mar-2022

